INSTRUCTIONS FOR USE

en

HM SPLINT

Hybrid Mesh Orthopedic Splint

Description:

The HM Splint is a synthetic splinting system, made from polyester fabric impregnated with a polyurethane resin and covered with span polar fleece on both sides. Exposure of HM Splint product to water initiates a chemical reaction which causes the product to become rigid.

Due to its flexibility, it easily achieves an elliptical form for articulated parts of the body, like an elbow, ankle, or knee, preventing injury from friction.

X-ray easily permeates and helps verify the detailed healing process of bones.

The finished HM Splint product is water repellent, it can be dried after showering with good ventilation and body heat.

Composition:

The HM Splint materials are composed of polyester fabric impregnated with a water-activated polyurethane resin, covered with span polar fleece on both sides.

Indications:

The HM Splint pre-cut and roll type products are intended for rigid immobilization through splinting procedures for the body's extremities, the immobilization of bone fractures, and the treatment of soft tissue and joint injuries.

Specific application suitability should be the responsibility of a qualified, on-site medical professional.

Specifications:

Package contents: 1 HM Splint

HM Splint is available in various sizes for use on the body's extremities.

Product type	Product Code	Size/ Variation/Indication			
1	HM SPLINT	HMS330	3 inch (7,62 cm x 30 cm)	Pre-cut type	One Step (Kid/ Adult)
2	HM SPLINT	HMS388	3 inch (7,62 cm x 88 cm)	Pre-cut type	
3	HM SPLINT	HMS438	4 inch (10,16 cm x 38 cm)	Pre-cut type	
4	HM SPLINT	HMS476	4 inch (10,16 cm x 76 cm)	Pre-cut type	
5	HM SPLINT	HMS576	5 inch (12,7 cm x 76 cm)	Pre-cut type	
6	HM SPLINT	HMS5114	5 inch (12,7 cm x 114 cm)	Pre-cut type	
7	HM SPLINT	HMS3450	3 inch (7,62 cm x 450 cm)	Roll type	
8	HM SPLINT	HMS4450	4 inch (10,16 cm x 450 cm)		
9	HM SPLINT	HMS5450	5 inch (12,7 cm x 450 cm)		

HM SPLINT Application procedure:


■ Short Arm Splint

1. If using a roll type, open the box, withdraw the plastic seal, then select and cut the most proper amount of the HM Splint. Place back the plastic seal to conserve the rest of the product.
2. If using a pre-cut type, select the most proper size for the patient. Remove splint from the pouch.
3. Submerge the pad into the water (20°C to 25°C) for 7 seconds.
4. Gently eliminate its water excess using a towel, then place the HM Splint pad on the part of the body to be treated. Fix it using a bandage.
5. Place the patient in the required position and maintain it still for 10 minutes.

Instruction for use:

Visit www.hmcast.eu for more information on the procedure.

MD Medical device

PRECAUTIONS AND NOTES:

Avoid breaking, folding, or perforating the packing before its use.
 Do not squeeze or wring. Do not use it if the packaging is damaged.
 Do not apply if show allergies, anaphylaxis, bleeding, inflammation, eruption, or infection.
 When using the HM Splint, it is necessary to wear protective gloves.
 Skin contact with uncovered polymer must be avoided. If the resin does come into contact with the skin, it can be removed with alcohol or acetone.

Shelf Life and Storage Requirements:

Store the product in a cool and dry place, away from the sunlight.
 Store it horizontally. If the product is stored in a standing position, the hardening liquid inside the product is withdrawn and the strength of the product is uneven, making it unusable.
 Keep at temperatures between 15°C and 30°C (ideally 22°C to 25°C.)
 For Shelf-life, refer to the expiry dates. Each box and pouch is marked with lots and expiry dates.

Disposal:

HM Splint can be incinerated or landfilled following local regulations.

Make sure there is good ventilation:

The coated polyester yarn with polyurethane resin contains a low volatile form of isocyanates.
 It is generally recommended to process the HM Splint in ventilated rooms.

Pay attention to the water temperature of 21°C -24°C when using:

For processing the HM Splint, clean water with a temperature of 21°C -24°C should be used.
 The setting reaction is exothermic, instruct the patient that the splint might feel warm during setting.
 This subsides within a few minutes. If warmer water (> 24°C) is used, the exothermic reaction is stronger; it can lead to a higher temperature development and undesired side effects. When processing the HM Splint, it is generally recommended to maintain the water temperature as specified in these instructions.

Explanation of the symbols on the label:

	Manufacturer		Keep away from sunlight
	Date of Manufacture		Upper and lower temperature limitation
	Use by date		Do no reuse
	Batch code		Consult instructions for use
	Catalog number		Caution
	Do not use if package is damaged		Keep dry
	Quantity in package		Authorized EC Representative in the European Community
	Product is not made with natural rubber		Medical Device

New Cast Industry Co., Ltd,
 11, Cheoyongsaneop 5-gil,
 Cheongnyang-myeon,
 Ulju-gun, Ulsan,
 Korea

Sanusmed Distribution Ltd,
 Calea Floreasca, no.60, 3rd floor,
 014462, Bucharest, Romania
 Phone: +40 721 261 363 / Fax: +40 212311031
 Email: office@hmcast.eu / Web: www.hmcast.eu

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